

PRODUCT DATA SHEET

PERFUSIA-SR®

DESCRIPTION

Dietary Supplement

Perfusia-SR®

Supplement Facts

Serving Size: Two Capsules Servings Per Container: 60

Two Capsules Contain:	% DV
L-Arginine **	1000 mg. *

*Daily Value (DV) not established

Other Ingredients: High and low viscosity hydroxypropyl methylcellulose and Magnesium Citrate/Laurate, Hypromellose (derived from cellulose) capsule, Leucine,[†] Silicon Dioxide.[†]

**Arginine complexed with Time-Sorb® - a patent-pending, time release matrix.

[†]A hypoallergenic ingredient used to fill space.

SA525 / 120 Vegetarian Capsules

EXCIPIENT FREE

This product does not contain wheat, gluten, corn, yeast, soy, egg, dairy products, or artificial colors, sweeteners, or flavors. This product also does not contain lactose, palmitic acid, or magnesium, calcium, or vegetable stearates.

HOW SUPPLIED

Perfusia-SR® is supplied in bottles of 120 two-piece vegetarian capsules. Packed 12 bottles per case.

CONTRAINDICATIONS AND WARNINGS

Individuals with known herpes infection should not take arginine without also taking lysine, as arginine alone has been noted to cause herpes eruptions in susceptible individuals.

Since polyamines act as growth factors for cancers, and arginine may stimulate polyamine synthesis, chronic administration of arginine in cancer patients should be avoided. L-arginine is contraindicated in hepatocellular carcinoma.

This product is contraindicated in a patient with a history of hypersensitivity to any of its ingredients.

INTERACTIONS

Perfusia-SR used concomitantly with Simvastatin promotes beneficial endothelium-dependent vasodilation in patients with elevated asymmetrical dimethylarginine (ADMA).

PREGNANCY WARNING

A patient who is pregnant or nursing, or who might become pregnant, should consult a health-care practitioner before using this product.

DOSAGE

Adult patients: 2 capsules twice daily, or as otherwise directed by a health-care practitioner.

STORAGE

Store tightly sealed in a cool dry place, away from direct light. Keep out of the reach of children.

TRADEMARK AND PATENT INFORMATION

Time-Sorb® is a patent-pending time released matrix.

MANUFACTURED BY

Thorne Research, Inc.
Dover, Idaho USA 83825
800-228-1966
208-265-2488 fax
info@thorne.com
www.thorne.com

*While some companies say
they are raising their standards
... at Thorne, we set them.™*



*the Original Specialists
in Pure Encapsulations™*

EMAIL: info@thorne.com WEBSITE: www.thorne.com
TOLL FREE: (800)228-1966 TOLL-FREE FAX: (800)747-1950
P.O. BOX 25, DOVER, IDAHO, 83825

These statements have not been evaluated by the Food and Drug Administration. This product is not intended to diagnose, treat, cure or prevent any disease.

© 2009 Thorne Research, Inc. All rights reserved.